


The software is easily accessed and viewed with the use of a web browser such as Internet Explorer or Mozilla Firefox. Once logged-in, the administrators have access to detailed graphs, charts, and reports, of performed transactions. The admin panel is also used to setup users and configure the server. Each server offers great flexibility and can be configured to fit particular pharmacies needs and requirements. For example, some servers may be configured in a very strict manner in order to track and control all users who perform transactions, whereas other servers may be configured more freely where users are not required to log in, and the server is mainly used for sharing drug data and information.



The screenshot displays the 'DRX-SSX SYSTEM' administrative interface. At the top, there is a navigation bar with tabs for 'Status', 'Reports', 'Users', 'Backup', 'Update', 'Config', and 'Database maintenance'. Below this, a sub-navigation bar highlights 'Scenes', 'Server', and 'Database'. The main content area is titled 'Networking and IP Configuration' and includes fields for 'Public port IP Address' (192.168.1.199) and 'Private port IP Address' (192.168.138.1), each with 'Automatic' and 'Manual' radio buttons. A 'Scenes listening port' field is set to '1000'. Below this is a 'Security Rules' section with checkboxes for 'Allow scales to connect on public ports', 'Allow anonymous login', 'Allow API Update by Anonymous users', and 'Allow transaction aborting'. The 'Transaction Reports' section has buttons to 'Delete all transactions older than 6 months' and 'Delete all transaction report logs'. The 'Report Tab View' section contains checkboxes for various report columns: 'Display LOT Column', 'Display Strength Column', 'Display User Column', 'Display Rx Size Column', 'Display Verification Column', 'Display Unit Name Column', and 'Display Unit ID Column'. A 'UDP Configuration' section is partially visible at the bottom.

< [Previous](#) | [Next](#) >